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

STEP-BY-STEP INSTRUCTIONS ON HOW TO PUT THE ORTHOSIS ON

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EN Before application

Open the green and white marked quick release of the brace. Push the quick release forward until you can lift it out through the keyhole (picture 1 and 2).

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

EN Application

Bend the leg to about 30°. Pull the brace up until the joint centre is approximately at the level of the upper edge of the kneecap (picture 3).

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

EN Make sure that the brace is straight and not twisted on the leg (picture 4).

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EN Now close the top, green marked quick release (picture 5) and secondly the low, white marked quick release (picture 6).

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

EN Make sure that the straps marked in white and green are attached to the corresponding attachment points on the hinge.

The Orthotist has preset the brace and explained the orientation value to you. If you require more or less tension adjust the length of straps by using the Velcro fasteners (picture 7).

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EN Take a few steps to verify the correct position on the leg. If necessary, optimize the position again (picture 8).

INSTRUCTION FOR USE VEOXGENU 1.0

Dear Customer,

thank you for choosing a veomedical product. Please read through these instructions for use carefully. If you have any questions, please contact your doctor or medical professional.

This instruction for use is also available as a PDF-Document and can be downloaded from www.veomedical.de. Visually impaired and blind people can also have the PDF read out to them using suitable software.

1. INTENDED PURPOSE

The veoXgenu 1.0 is a medical device. It's a brace for targeted unilateral relief and stabilisation of the knee joint.

2. INTENDED USER

This device is intended for people with unicompartmental knee joint osteoarthritis or meniscal tears.

3. INDICATIONS FOR USE

- › Mild to severe medial or lateral osteoarthritis
- › Degenerative meniscal tears
- › Cartilage bone lesions
- › Bone bruise
- › Other knee conditions that may benefit from unicompartmental unloading
- › Condition after cartilage surgery or cartilage build-up with necessary medial or lateral post-operative unloading.

4. CLINICAL BENEFIT

The brace provides a targeted relief and stabilization for more movement with osteoarthritis of the knee.

5. CONTRAINDICATIONS

Bilateral osteoarthritis (medial and lateral)

6. BODY SEGMENT

Knee joint

7. RISKS OF USING THIS PRODUCT

- › The product provides its unloading effect during physical activity and stretching of the leg.
- › For first time use, we recommend that you gradually increase the usage time while building up tolerance to the brace.
- › Remove your brace during prolonged periods of rest.
- › Use the brace only as indicated and in accordance with the further instructions of the medical professional.

- › If used with other products, consult your medical professional.
- › Do not make any accidental changes to the product, which may alter the effect and cause health damage. Warranty and liability are excluded in such cases.
- › Avoid contact with fatty and acidic products or lotions.
- › Outside applied braces may cause local pressure symptoms or, in rare cases, constrict passing blood vessels or nerves.
- › If you notice changes or increasing discomfort while wearing the product, stop further use and contact your medical professional.
- › In the following clinical pictures, the use of the product is only indicated after consulting your doctor:
 - Skin diseases/injuries in the treated part of the body, especially in cases of inflammation, scars with swelling, redness and overheating
 - varicose veins (varicosis)
 - Strong Baker Cyst.

8. NOTES ON REUSE

This product is intended and tested for use on a single patient. Use on more than one patient is excluded.

9. ASSEMBLY INSTRUCTIONS

The product must be fitted by an Orthotist or qualified medical professional.


- Adjustment of the flexion and extension angle (to be performed by the Orthotist only)
- Setting the extension and flexion limits
- Extension limiting wedge, 0°
- Flexion limiting wedge, 90°.

Assembly

An 0° extension wedge is preinstalled in your brace. To replace the wedges, loosen the screws on the inside and remove it.

Insert the extension wedge from the front and the appropriate flexion wedge from the back into the side opening and screw the wedge tightly.

10. APPLICATION INSTRUCTIONS

Measurement Guide

Measure the circumference of the thigh at 15 cm (6 inches) above the mid-patella and measure the circumference of the calf at 15 cm (6 inches) below the mid-patella with knee slightly bent.

➤ See page 6 to 10 for step-by-step instructions on how to put the orthosis on.

11. CLEANING INSTRUCTIONS

Insert the entire brace in the laundry bag without removing the straps.

- › Close Velcro fasteners before washing.
- › Hand wash at 30° with a gentle detergent.
Higher temperatures can destroy the glued seams.
- › Do not wash with softener
- › Do not bleach
- › Do not iron
- › Do not dry clean
- › Air dry, do not dry in the dryer.

i NOTES

Avoid contact with chlorinated or salt water.
In case of contact rinse thoroughly with tap water.

12. PRODUCT MATERIAL CONTENT

Hinge: Aluminium, POM

Textile sleeve: 59% Polyamid, 41% Elastane

Strap: 100% Polyamide

Pads: ABS

Adhesive edge: 100% Silicone

Adhesive materials: 100% Polyurethane

13. STORAGE INSTRUCTIONS

Please store the product dry and protect it from direct sunlight.

14. WARRANTY

The legal provisions of the country in which you purchased the product apply. If you suspect a warranty case, please first contact your medical professional directly from whom you received the brace. Please note the instructions for use and care, otherwise the warranty may be impaired

or excluded. The warranty is excluded in the event of use that is not in accordance with the indications, in the event of non-compliance with the instructions of the medical professional and in unauthorized changes to the product.

15. DUTY TO REPORT

Due to legal requirements, you are obliged to report any serious incident when using this brace both to the manufacturer and to the BfArM (Federal Institute for Drugs and Medical Devices) immediately.

16. ENVIRONMENT/DISPOSAL

All product and packaging components should be disposed of in accordance with the applicable environmental regulations of the respective country. If you have any questions about recycling or environmentally sound disposal, please contact the competent authorities.

17. MANUFACTURER

veomedical GmbH

Enzianweg 8

50259 Pulheim

e-mail: hallo@veomedical.de

www.veomedical.de

Hersteller

Manufacturer:

veomedical GmbH

hallo@veomedical.de

www.veomedical.de

veomedical.de