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EN Before application

Open both straps (picture 1) and put the cooling compress to the inner side of the brace.

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
EN Application

Place the orthosis on the leg from the front so that the lower edge of the kneecap recess sits on the lower edge of the kneecap (picture 2).

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
EN First wrap the lower strap around the leg and close it (picture 3). Check the fit of the orthosis.

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
EN Then wrap the upper strap around your leg and close it too. Check that the orthosis fits well anywhere on the leg and does not slip (picture 4).

If the orthosis slips, close the straps a little tighter.

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
EN Brace removal

Open the two straps. Pull the brace upwards (picture 5).

5

INSTRUCTION FOR USE

OPPO GENU COOL

Dear Customer,

thank you for choosing an OPPO Medical product. Please read through these instructions for use carefully. If you have any questions, please contact your doctor or medical professional.

This instruction for use is also available as a PDF-Docum-ent and can be downloaded from www.oppomedical.de. Visually impaired and blind people can also have the PDF-Docum-ent read out to them using suitable software.

1. INTENDED PURPOSE

The OPPO GENU COOL is a medical device. It is an orthosis for stabilizing and guiding the knee joint with an optional cooling function.

2. INTENDED USER

This device is intended for people with indications that require reinforcement of the physiological stabilisation of the knee joint.

3. INDICATIONS FOR USE

All indications that require reinforcement of the physio-logical stabilisation of the knee joint and/or stress relief of the joint, such as:

- › injury to the collateral ligament
- › meniscal injury
- › surgical intervention
- › knee instability

4. CLINICAL BENEFIT

The brace provides reinforcement of the physiological stabilisation and guidance of the knee joint.

5. CONTRAINDICATIONS

None known at present.

6. BODY SEGMENT

Knee joint.

7. PERFORMANCE CHARACTERISTICS

The brace adopts a series of advanced product designs. The polycentric joint, together with the straps and the textile, create a stabilizing and guiding effect on the knee joint.

8. WARNING / PRECAUTION

- › Use the brace only as indicated and in accordance with the further instructions of the medical profes-sional.

- › If used with other products, consult your medical professional.
- › Do not make any accidental changes to the product, which may alter the effect and cause health damage. Warranty and liability are excluded in such cases.
- › Avoid contact with fatty and acidic products, ointments or lotions.
- › Outside applied braces may cause local pressure symptoms or, in rare cases, constrict passing blood vessels or nerves.
- › If you notice changes or increasing discomfort while wearing the product, stop further use and contact your medical professional.
- › Regular practitioner supervision is recommended for users with:
 - Skin diseases/injuries in the treated part of the body, especially in cases of inflammation, scars with swelling, redness and overheating
 - Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
 - Lymph drainage disorders – equally
 - swelling of soft tissue adjacent to the area of application which is not related to a condition
- › Do not use this device if there is any damage on the device itself.

9. NOTICE TO ADDITIONAL COOLING FUNCTION

Depending on the doctor's instructions, the orthosis can also be used with the included cooling compress.

To use the cooling function, leave the included cooling compress in the freezer briefly (max. 30 minutes) or simply in the fridge. An excessively long cooling period in the fridge is not recommended as the temperature could then drop too low and frostbite could occur on application.

The cooling compress can simply be fastened to the black inner side of the brace. The blue inner side of the cooling compress is left facing you.

10. WARNING / PRECAUTION TO THE COOLING FUNCTION

Please consult the doctor treating you before applying cold therapy. This is mandatory. It may be that this treatment is contraindicated for you.

Ensure that an underlay (thin towel or piece of clothing) is present between the skin and the cooling compress. The length of time and frequency of the cooling phases must always be specified by the treating doctor and should not last longer than 30 minutes.

If the cooling compress is defective for any reason, if it leaks and the gel comes into contact with the skin or eyes, wash the affected areas thoroughly with water straight away. Never try to heat the cooling compress in any way (e.g. water baths, microwaving, ovens etc. are all prohibited).

11. CONTRAINDICATIONS FOR COOLING THERAPY

If you have a vascular disease, circulation disorder, paraesthesia or tissue inflammation, these therapies may only be introduced once they have been discussed with your doctor.

12. NOTES ON REUSE

This product is intended and tested for use on a single patient. Use on more than one patient is excluded.

13. ASSEMBLY INSTRUCTIONS

The product must be fitted by an Orthotist or qualified medical professional.

Position both of the joints on the brace. To do this, put the joint splints in the small attached hook-and-loop pockets.

Now stick the splints onto the side of the brace. Ensure that the brace joints are in the correct position:

The middle of the brace joints should be level with the top of the kneecap

The joints must be positioned behind the centre of the side of the leg and must run parallel to each other. The joints must be adapted to the patient's individual leg shape by limiting them.

If the straps are too short, simply fasten both of the additional strap extensions on the hook-and-loop areas at the ends of the straps for the brace. You can now wrap the straps around the leg as usual and fasten them on the brace.

14. APPLICATION INSTRUCTIONS

One size fits all. Maximum circumference on the thigh: 65 cm without belt extension – 90 cm with belt extension.

➤ See page 6 to 9 for step-by-step instructions on how to put the orthosis on.

12. CLEANING INSTRUCTIONS

- › Close Velcro fasteners before washing
- › Hand wash at 30° with a gentle detergent
- › Do not wash with softener
- › Do not bleach
- › Do not iron
- › Do not dry clean
- › Air dry, do not dry in the dryer

i NOTES

Avoid contact with chlorinated or salt water. In case of contact rinse thoroughly with tap water.

16. PRODUCT MATERIAL CONTENT

Knee sleeve: cotton, nylon

Cooling compress: propylene glycol, sodium carboxymethyl cellulose, water

Joints: aluminium, POM

17. STORAGE INSTRUCTIONS

Please store at room temperature. Avoid sunlight and/or exposure to excessive heat or moisture.

18. WARRANTY

The legal provisions of the country in which you purchased the product apply. If you suspect a warranty case, please first contact your medical professional directly from whom you received the brace. Please note the instructions for use and care, otherwise the warranty may be impaired or excluded. The warranty is excluded in the event of non-indication-compliant use, non-compliance with the instructions of the medical professional and in the event of unauthorized product modification.

19. DUTY TO REPORT

Serious incidents that occur in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user resides. Our contact details can be found on this brochure.

20. ENVIRONMENT / DISPOSAL

All product and packaging components should be disposed of in accordance with the applicable environmental regulations of the respective country. If you have any questions about recycling or environmentally sound disposal, please contact the competent authorities.

Hersteller Manufacturer:

OPPO MEDICAL INC.
Seattle City Center
1420 Fith Ave., Ste. 220080
Seattle WA 98101, USA
info@oppomedical.com
www.oppomedical.com
www.oppomedical.de

**EU Bevollmächtigter
EU Authorized representative:**

MT Promedt Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert, Germany
info@mt-procons.com
www.mt-procons.com

EU Importeur EU Importer:

veomedical GmbH
hallo@veomedical.de
www.veomedical.de

Artikelnummer Article Number:
324D100655U0

Datum Date:
09.2023, OR Version