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

EN Before application

Open the two front velcro fasteners of the brace without removing the straps from the slots of the joint splints (Picture 1).

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EN Before Application


Open the two front velcro fasteners. Close the hook straps loosely in the middle of the fleece straps. This means you can put the orthosis on comfortably without causing damage to the 3D fabric (Picture 2).

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
EN Application

Pull the brace over the affected knee until the center of the joints are approximately at the middle of the kneecap (Picture 3). If you have limited grip strength, use the two inner finger loops for correct positioning (Picture 4).

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EN Make sure that the brace is not twisted on the leg. If the green patterned flex area is directly over the kneecap, the orthosis is correctly positioned (Picture 5).

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
EN Tighten the two front straps until you achieve a comfortable level of stabilization and close the Velcro fasteners (Picture 6).

The back straps have already been adjusted by your medical professional. If you need to readjust them, proceed as you did with the front straps.

INSTRUCTION FOR USE

OPPO GENU S

Dear Customer,

thank you for choosing an OPPO Medical product. Please read through these instructions for use carefully. If you have any questions, please contact your doctor or medical professional.

This instruction for use is also available as a PDF-Docum-ent and can be downloaded from www.oppomedical.de or veomedical.de. Visually impaired and blind people can also have the PDF-Docum-ent read out to them using suitable software.

1. INTENDED PURPOSE

The OPPO GENU S is a medical device. It is an orthosis for stabilizing and guiding the knee joint.

2. INTENDED USER

This device is intended for people with indications that require reinforcement of the physiological stabilisation of the knee joint.

3. INDICATIONS FOR USE

All indications that require reinforcement of the physio-logical stabilisation and/or guidance of the knee joint, such as:

- injury to the anterior cruciate ligament, lateral collateral ligament, medial collateral ligament
- meniscal injury
- surgical intervention
- knee instability.

4. CLINICAL BENEFIT

The brace provides reinforcement of the physiological stabilisation and guidance of the knee joint.

5. CONTRAINDICATIONS

None known at present.

6. BODY SEGMENT

Knee joint.

7. PERFORMANCE CHARACTERISTICS

The brace adopts a series of advanced product designs. The anatomical joints, together with the straps and the 3D knit textile, create a stabilizing and guiding effect on the knee joint.

8. WARNING / PRECAUTION

- › Use the brace only as indicated and in accordance with the further instructions of the medical professional.

- › If used with other products, consult your medical professional.
- › Do not make any accidental changes to the product, which may alter the effect and cause health damage. Warranty and liability are excluded in such cases.
- › Avoid contact with fatty and acidic substances, creams, or lotions.
- › Outside applied braces may cause local pressure symptoms or, in rare cases, constrict passing blood vessels or nerves.
- › If you notice changes or increasing discomfort while wearing the product, stop further use and contact your medical professional.
- › Regular practitioner supervision is recommended for users with:
 - skin diseases/injuries in the treated part of the body, especially in cases of inflammation, scars with swelling, redness and overheating
 - paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
 - lymph drainage disorders
 - swelling of soft tissue adjacent to the area of application which is not related to a condition.
- › Do not use this device if there is any damage on the device itself.


9. NOTES ON REUSE

This product is intended and tested for use on a single patient. Use on more than one patient is excluded.

10. ASSEMBLY INSTRUCTIONS

The product must be fitted by an Orthotist or qualified medical professional. When selecting, the correct side and the correct size must be determined. The straps can be cut to the appropriate length with scissors.

11. APPLICATION INSTRUCTIONS

Measurement Guide

Measure the circumference of the thigh at 15 cm (6 inches) above the mid-patella and measure the circumference of the calf at 15 cm (6 inches) below the mid-patella with knee fully extended.

Adjusting the brace

Open all four straps and position the brace. Adjust the lengths of the four straps, ensuring that the two joint splints follow the leg anatomy even after adjusting the straps. Make sure that the brace is not twisted. The straps can be shortened at the back end with scissors.

- See page 5 to 7 for step-by-step instructions on how to put the orthosis on.

12. CLEANING INSTRUCTIONS

- › Close velcro fasteners before washing
- › Hand wash at min. 30°C with a gentle detergent
- › Do not wash with softener
- › Do not bleach
- › Do not iron
- › Do not dry clean
- › Air dry, do not dry in the dryer.

NOTES

Avoid contact with chlorinated or salt water. In case of contact rinse thoroughly with tap water.

13. PRODUCT MATERIAL CONTENT

Knee sleeve: 76% Polyamide, 13% Polyester,
11% Elastane

Straps: 100% Polyester

Joints: 100% Polyamide

14. STORAGE INSTRUCTIONS

Please store at room temperature. Avoid sunlight and/or exposure to excessive heat or moisture.

15. WARRANTY

The legal provisions of the country in which you purchased the product apply. If you suspect a warranty case, please first contact your medical professional directly from whom you received the brace. Please note the instructions for use and care, otherwise the warranty may be impaired or excluded. The warranty is excluded in the event of non-indication-compliant use, non-compliance with the instructions of the medical professional and in the event of unauthorized product modification.

16. DUTY TO REPORT

Serious incidents that occur in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user resides. Our contact details can be found on this brochure.

17. ENVIRONMENT / DISPOSAL

All product and packaging components should be disposed of in accordance with the applicable environmental regulations of the respective country. If you have any questions about recycling or environmentally sound disposal, please contact the competent authorities.

NOTIZEN
NOTES

Hersteller Manufacturer:

OPPO MEDICAL INC.
Seattle City Center
1420 Fith Ave., Ste. 220080
Seattle WA 98101, USA
info@oppomedical.com
www.oppomedical.com
www.oppomedical.de

**EU Bevollmächtigter
EU Authorized representative:**

MT Promedt Consulting GmbH
Altenhofstr. 80
66386 St. Ingbert, Germany
info@mt-procons.com
www.mt-procons.com

EU Importeur EU Importer:

veomedical GmbH
hallo@veomedical.de
www.veomedical.de

Artikelnummer Article Number:
324D101655U0

Datum Date:
09.2024, 0R Version