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
STEP-BY-STEP INSTRUCTIONS ON HOW TO PUT THE ORTHOSIS ON

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
EN Before application

Open the two front hook-and-loop fasteners of the orthosis without removing the straps from the slots of the joint bars. (Image 1)

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EN Before Application

Fasten the hook straps loosely, approximately in the middle of the loop straps. You can put on the orthosis comfortably without causing any damage to the 3D sleeve. (Image 2)

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
EN Application


Pull the orthosis over the affected knee until the center of the side joints is approximately at the same height as the middle of the kneecap. (Image 3)

If you have limited grip strength, use the two inner loops. (Image 4)

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EN Make sure the orthosis is not twisted on your leg. When the green patterned flex zone is positioned directly above the kneecap, the orthosis is correctly aligned. (Image 5)

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
EN Tighten the two front straps until you achieve a comfortable level of stabilization and close the hook-and-loop fasteners. (Image 6)

The back straps have already been adjusted by your qualified orthopaedic technician. If you need to readjust them, proceed in the same way as with the front straps.

INSTRUCTION FOR USE

OPPO GENU S+

Dear Customer,

thank you for choosing an OPPO Medical product. Please read through these instructions for use carefully. If you have any questions, please contact your doctor or medical professional.

These instructions for use are also available in electronic form as an accessible PDF document and can be downloaded from www.veomedical.de. The document is designed to be compatible with assistive technologies (e.g. screen readers). A version of the instructions for use in easy-to-understand language is available upon request.

1. INTENDED PURPOSE

The OPPO GENU S+ is a medical device. It is a knee orthosis with lateral joint bars designed to stabilize and support the knee joint, allowing limitation of flexion and extension as required for therapy.

2. INTENDED USER

This medical device is intended for individuals with mild to moderate knee complaints requiring stabilization and/or guidance of the knee joint.

3. INDICATIONS FOR USE

All indications that require reinforcement of the physiological stabilisation of the knee joint and/or stress relief of the joint, such as:

- › injury to the collateral ligament
- › meniscal injury
- › surgical intervention
- › knee instability (e.g. knee osteoarthritis)

4. CLINICAL BENEFIT

This medical device stabilizes and guides the knee joint and provides the option to restrict knee flexion according to medical necessity.

5. CONTRAINDICATIONS

None known at presence.

6. BODY SEGMENT

Knee joint.

7. PERFORMANCE CHARACTERISTICS

The orthosis utilizes several advanced product features. The anatomically shaped joint bars, in combination with circumferential straps and 3D-knitted fabric, provide a stabilizing and guiding effect on the knee joint. The range of motion in flexion and extension can be limited in line with therapy requirements.

8. WARNING / PRECAUTION


- › Use the orthosis only for the intended indication and in accordance with the advice of medical professionals.
- › When using other products simultaneously, consult your physician.
- › Do not modify the product on your own, as this may alter its effectiveness and cause injury. Warranty and liability are void in such cases.
- › Avoid contact with greasy or acidic substances, ointments or lotions.
- › Externally applied aids that are too tight may cause local pressure points or, in rare cases, compress blood vessels or nerves.
- › If you notice changes or increasing discomfort while wearing the product, discontinue use and consult your doctor.
- › In the following conditions, use of the product should only occur after consultation with a physician:
 - Skin disorders or injuries in the treated area, especially inflammatory, oozing or open lesions, or raised scars with swelling, redness or warmth.
 - Sensory or circulatory disorders (e.g. diabetes, varicose veins)
 - Lymphatic disorders
 - Unclear soft tissue swellings, even outside the treated area.
- › Do not use the orthosis if the product itself is damaged.

9. NOTES ON REUSE

This product is designed and tested for use by a single patient only. Reuse on another person is not permitted.

10. ASSEMBLY INSTRUCTIONS

The product must only be fitted by qualified orthopaedic technicians.

Determine the correct side and size.

Measure the thigh circumference 15 cm (6 inch) above and the calf circumference 15 cm (6 inch) below the center of the knee with the leg extended.

Adjusting the brace

Open all straps and position the orthosis. Adjust the strap lengths so that the joint bars follow the leg's natural contour and the brace is not twisted. The straps can be shortened at the rear ends with scissors.

Joint adjustment

Do not adjust the joints while the orthosis is being worn. The package includes extension wedges at 0° and 30°, and flexion wedges at 0°, 30°, 45°, 90°. Additional wedges are available upon request. All flexion and extension wedges can be inserted and removed without tools.

- › Select the required range limitation according to the physician's prescription. All wedges are marked with their respective angle.
- › Insert the flexion wedge on the dorsal side of the joint first, followed by the extension wedge on the ventral side. The wedges are correctly engaged when a distinct 'click' is heard.
- › Always use the same wedge limitations in both joints to ensure symmetrical restriction.
- › Check the secure fit of the wedges before applying the orthosis to the patient. To remove the wedges, pull firmly on the wedge tab. To avoid pinching of skin or clothing in the joint, always insert two wedges per joint – one flexion and one extension wedge.

Functional check

The orthosis must not slip or cause constriction. After 10–15 minutes of wear, perform a brief fit check.

11. APPLICATION INSTRUCTION

Wearing time should follow the medical professional's recommendation and the patient's condition. Initially, the orthosis should be worn for shorter periods (e.g. 1–2 hours per day) to allow the body to adapt. The duration can then be gradually increased as advised.

The expected service life of the product is approximately 6 months with proper use and care. After this period, the orthosis should be inspected by a professional for wear or functional changes and replaced if necessary.

Do not make unauthorized modifications to the joint mechanism, straps or textile components. Improper alterations may compromise effectiveness or cause injury.

- See page 6 to 8 for step-by-step instructions on how to put the orthosis on.

12. CLEANING AND DESINFECTION

Visible soiling must be removed prior to any disinfection process.

Cleaning

The product can be cleaned by hand washing using a mild detergent at a temperature of at least 30°C. Close all hook-and-loop fasteners before washing. Do not bleach, do not iron, and do not tumble dry.

Washer disinfection (preferred)

Cleaning may be performed using a disinfecting detergent or laundry sanitiser at a minimum temperature of 30°C. Only products suitable for textile materials shall be used.

Wipe or spray disinfection (alternative)

If washer disinfection is not possible, wipe or spray disinfection may be performed. Use a suitable, material-compatible surface disinfectant (e.g. Verbund für Angewandte Hygiene (VAH) -listed or equivalent). Ensure that all surfaces are thoroughly wetted. The contact time specified by the manufacturer of the disinfectant must be observed. Afterwards, allow the product to dry completely.

NOTES

- › Avoid contact with chlorinated or salt water. If contact occurs, rinse the orthosis thoroughly with tap water.
- › Only use material-compatible disinfectants
- › Do not use aggressive chemical agents
- › Do not sterilise the product
- › All components can be cleaned at temperatures of at least 30°C.

13. PRODUCT MATERIAL CONTENT

Knee sleeve:	76% polyamide, 13% polyester, 11% elastane
Straps:	100% polyester
Joint bars:	100% polyamide

14. STORAGE INSTRUCTIONS

Store the product at room temperature. Protect it from direct sunlight and moisture.

Hersteller Manufacturer:

OPPO MEDICAL INC.
Seattle City Center
1420 Fifth Ave., Ste. 220080
Seattle WA 98101, USA
info@oppomedical.com
www.oppomedical.com
www.oppomedical.de

**EU Bevollmächtigter
EU Authorized representative:**

MT Promedt Consulting GmbH
Ernst-Heckel-Str. 7
66386 St. Ingbert, Germany
info@mt-procons.com
www.mt-procons.com

EU Importeur EU Importer:

veomedical GmbH
Enzianweg 8
50259 Pulheim, Germany
hallo@veomedical.de
www.veomedical.de

Artikelnummer Article Number:
324D102655U0

Datum Date:
03.2026, 1R Version

oppomedical.com